

For Research Use Only.
Not For Use In Diagnostics.

CoraLite594-conjugated Mouse Anti-Heavy Chain of Rabbit IgG

Catalog Number:SA00014-5



www.ptglab.com

Information

Catalog Number:

SA00014-5

Host:

Mouse

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

CoraLite®594

Recommended Dilutions

1:100-1:500

Preservative

Proclin300

Fluorophore

CL594 Fluorescent Dye. Amax=593 nm, Emax=614 nm

Stabilizer

0.5% BSA

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

Purity

Protein G purification

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data