

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite647-conjugated AffiniPure F(ab')<sub>2</sub> Fragment Goat Anti-Rabbit IgG (H+L)



Catalog Number: SA00014-9

## Information

Catalog Number:

SA00014-9

Host:

Goat

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

CoraLite®647

## Recommended Dilutions

1:100 – 1:500

## Preservative

0.01% Sodium Azide

## Fluorophore

CL647 Fluorescent Dye. Amax=642nm, Emax=662nm

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A651/A280 = 4.2

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

