

# Rabbit anti-Human IgG

Catalog Number: B900110

---

Physical State	Liquid
Protein Concentration	1 mg/ml
Size	100ug
Purification	The antibody was purified from antisera by immunofinity chromatography.
Buffer	PBS, 50% glycerol, pH 7.3
Preservative	0.02% Sodium Azide
Storage	-20°C, avoid repeated thawing and freezing
Note	This product is for research use only, not for diagnostic or therapeutic use.

---

For technical support and original validation data for this product please contact:

T: 027-87531629    E: Proteintech-CN@ptglab.com    W: ptgcn.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.