

For Research Use Only.



# Goat anti-Human IgG

Catalog Number: B900220

---

<b>Physical State</b>	Liquid
<b>Protein Concentration</b>	1 mg/ml
<b>Size</b>	100ug
<b>Purification</b>	The antibody was purified from antisera by immunoaffinity chromatography.
<b>Buffer</b>	PBS, 50% glycerol, pH 7.3
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	-20°C, avoid repeated thawing and freezing
<b>Note</b>	This product is for research use only, not for diagnostic or therapeutic use.

---

For technical support and original validation data for this product please contact:

T: 027-87531629 E: Proteintech-CN@ptglab.com W: ptgcn.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.