

For Research Use Only

anti-p53 (N-term) recombinant VHH, unconjugated



www.ptgcn.com

Catalog Number: pt

Basic Information

Catalog Number:

pt

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2631380

Molecular Weight:

12.96 kDa

Description

Alpaca anti-p53 VHH, purified recombinant binding protein

Specificity/Target

p53

Species-Reactivity: Tested on human

Isoform: Alpha, beta and gamma

Epitope: 1-81 aa

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Structure of VHH.