

Catalog Number: **s3tma**

Basic Information

Catalog Number:
s3tma

Applications:
IP, Co-IP

Conjugate:
Magnetic Agarose beads; ~40 um (cross-linked 6% magnetic agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant

Description

The ChromoTek SUMO3-Trap Agarose consists of an anti-SUMO3 VHH, which is coupled to agarose beads. It can be used for the immunoprecipitation of SUMO3 and SUMO2 proteins from cell extracts of various organisms.

Specificity/Target

Binds specifically to human SUMO2 and SUMO3 (Small Ubiquitin-like Modifier 2 and 3). The ChromoTek SUMO3-Trap Agarose will immunoprecipitate both monovalent SUMO2/3 and chains formed by the process of SUMOylation. In addition to human SUMO2 und SUMO3, it will also bind to homologs from other mammals.

Elution buffer

2x SDS-sample buffer (Lämmli), 200 mM glycine pH 2.5

Affinity (K_D)


/

Storage

Storage:
Shipped at ambient temperature. Upon receipt store at +4°C. Stable for one year. Do not freeze!

Storage Buffer:
20% ethanol

Selected Validation Data



The SUMO3-Trap Magnetic Agarose (s3tma) was used to immunoprecipitate endogenous SUMO3ylated and SUMO2ylated proteins from a human cell line. Samples of the input lysate (I), non-bound flow-through (F) and bound (B) fraction were analysed using Western blot. SUMO2/3ylated proteins are highly enriched using the SUMO3-Trap. Data kindly provided by Prof. Alfred Vertegaal, Leiden University, The Netherlands.