

For Research Use Only

TurboGFP RecAb, VHH-rabbit IgG Fc fusion



www.ptgcn.com

Catalog Number: **tbfrb**

Basic Information

Catalog Number:
tbfrb

Clone No.:
CTK0204

Applications:
IF

Host/IsoType:
Rabbit/IgG

Conjugate:
Unconjugated

Type:

Primary antibody; chimeric fusion of Nanobody (VHH) and rabbit IgG Fc domain

Class:

Recombinant

RRID:

AB_2889371

Purification Method:

Recombinantly expressed and purified by IMAC

Description

Rabbit monoclonal antibody to TurboGFP. Recombinant fusion of anti-TurboGFP VHH [CTK0204] and rabbit IgG Fc domain.

Specificity/Target

TurboGFP

Recommended Dilution

IF/ICC: 1:400 - 1:2,000

Storage

Storage:

Shipped on ice. Aliquot upon receipt and store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:

PBS, 10% glycerol, Preservative: 0.02 % sodium azide, safety datasheet (SDS): sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells transiently expressing Actin-Chromobody2-TurboGFP were immunostained with TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204] (tbfrb-20, 1:1,000) and Nano-Secondary? alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor? 647 [CTK0101, CTK0102] (srbAF647-1, 1:1,000). Scale bar, 10 μ M.