

Catalog Number: v5ab

5 Publications

## Basic Information

**Catalog Number:**  
v5ab

**Clone No.:**  
SV5-P-K

**Applications:**  
WB

**Host/IsoType:**  
Mouse/IgG1 kappa

**Conjugate:**  
Unconjugated

**Type:**  
Primary Antibody

**Class:**  
Monoclonal

**RRID:**  
AB\_2868500

**Purification Method:**  
Protein A purified

## Description

Mouse monoclonal [SV5-P-K] to V5-tag

## Specificity/Target

V5-Tag (GKPIPNPLLGLDST)

## Recommended Dilution

Western Blot: 1:1,000

## Storage


**Storage:**

Shipped on blue ice. Store at -20°C/-4°F. Stable for 1 year after shipment.


**Storage Buffer:**

PBS, preservative: 0.02 % Proclin

Selected Validation Data



Western Blot (WB) analysis of His-eGFP-V5 added to HEK293T cell extract



The anti-V5 antibody [SV5-P-K] is a recombinant, monoclonal mouse IgG1, has very low background and detects even low amounts of V5-tagged proteins. Secondary antibody ChromoTek's Alpaca anti-mouse IgG1 Alexa Fluor 488 Nano-Secondary.