

For Research Use Only

anti-Vimentin recombinant VHH, unconjugated



www.ptgcn.com

Catalog Number: vt

Basic Information

Catalog Number:

vt

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type:

Nanobody

Class:

Recombinant

RRID:

AB_2827563

Molecular Weight:

26.7 kDa

Description

Alpaca anti-Vimentin VHH, purified recombinant binding protein

Specificity/Target

Vimentin

Vimentin-Label, a Vimentin VHH conjugated to ATTO488, specifically binds to the vimentin intermediate filament protein. Validated on canine MDCK cells, rodent BHK cells and on human A549 cells upon stimulation with TGF β .

Storage

Storage:

Shipped on blue ice. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Structure of Vimentin.